


EPR-spectroscopy for detection of free radicals and reactive oxygen species




Free radicals are particles possessing unpaired electrons, hence they are paramagnetic and highly reactive. They can be neutral or charged (ion-radicals), with one or more unpaired electrons (polyradicals), short-living (fractions of a second) or long-living (up to several years lifetime) at 298 K, solid, liquid or gaseous.

Free radicals play a vital role in redox, photochemical [1-3] and catalytic [4, 5] processes, as well as in main industrial processes: polymerization, telomerization, pyrolysis, cracking, heterogeneous catalysis. In order to stabilize the radicals, low temperature techniques are widely used, for example – liquid nitrogen cooling.

Long-living (stable) free radicals are paramagnetic compounds with highly delocalized free electrons and sterically hindered reaction centers [6].



Spin trapping is one of the most important methods in biological and medical research since it allows to detect and identify free radicals generated in living cells and tissues. Among them are the superoxide and hydroxyl radicals (the so-called Reactive Oxygen Species, ROS) as well as nitrogen oxide [7].



Free radicals and particularly – ROS formation and degradation are important chemical processes occurring in human body. Spin trapping technique allows to study the antioxidant properties of various substances and to detect ROS as well as some other types of free radicals [8]. Various chemical reagents (PBN, DMPO) are used as spin traps [9].

Figure 1 – Ascorbyl and DMPO-OH radical EPR spectrum (radical was generated during the Fenton-like reaction with Cu²⁺)





Figure 2 – DMPO-OH and DEPMPO-OH EPR spectra (radicals were generated by UV-irradiation)



REFERENCES

1. Z. Wang, V. Srivastava, S. Iftekhar, I. Ambat, M. Sillanpää, Fabrication of Sb₂O₃/PbO photocatalyst for the UV/PMS assisted degradation of carbamazepine from synthetic wastewater, *Chemical Engineering Journal* (2018), doi: <https://doi.org/10.1016/j.cej.2018.08.068>
2. Z. Wang, A. Deb, V. Srivastava, S. Iftekhar, I. Ambat, M. Sillanpää. - Investigation of textural properties and photocatalytic activity of PbO/TiO₂ and Sb₂O₃/TiO₂ towards the photocatalytic degradation Benzophenone-3 UV filter. - *Separation and Purification Technology*, 228 (2019) 115763-115775.
3. B. Gao, S. Iftekhar, V. Srivastava, B. Doshi, M. Sillanpää. - Insights into the generation of reactive oxygen species (ROS) over polythiophene/ZnIn₂S₄ based on different modification processing. *Catal. Sci. Technol.*, 2018, DOI: 10.1039/C8CY00303C.
4. Ren, Z., Romar, H., Varila, T., Xu, X., Wang, Z., Sillanpää, M., Leiviskä, T., Ibuprofen degradation using a Co-doped carbon matrix derived from peat as a peroxymonosulphate activator, *Environmental Research*, <https://doi.org/10.1016/j.envres.2020.110564>.
5. P. Kaur, S. Frindy, Y. Park, M. Sillanpää, M.A. Imteaz. - Synthesis of Graphene-Based Biopolymer TiO₂ Electrodes Using Pyrolytic Direct Deposition Method and its Catalytic Performance. - *Catalysts* 2020, 10, 1050-1066.
6. Rozantsev, E.G., Goldfein M.D., Pulin V.F. - Organic paramagnetics. - Saratov State University, 2000. - 340 p.
7. S. Suzen, H. Gurer-Ohran and L. Saso. - Detection of Reactive Oxygen and Nitrogen Species by Electron Paramagnetic Resonance (EPR) Technique *Molecules* 2017, 22, 181.
8. Q. Wang, D. Zheng, M.E. McKinnon, X.-Q. Yang and D. Qu - Kinetic Investigation of Catalytic Disproportionation of Superoxide Ions in the Nonaqueous Electrolyte used in Li-Air Batteries. - BNL-107735-2015-JA.
9. J.M. Fontmorin, R.C. Burgos Castillo, W.Z. Tang, M. Sillanpää. - Stability of 5,5-dimethyl-1-pyrroline-N-oxide as a spin-trap for quantification of hydroxyl radicals in processes based on Fenton reaction. - *Water Research* 99 (2016) 24-32.

Some examples of using benchtop EPR spectrometer SpinScan X for free radicals detection:

- 1) Ascorbyl radical is generated during the ascorbic acid oxidation – activation by virtually any reactive species occurring in biological media (fig.1a).
- 2) hydroxyl radicals detection during the Fenton-like reaction with Cu²⁺. Generation at pH 6.0, 1 mM Cu²⁺, 1 mM ascorbic acid and 1 mM H₂O₂ with 100 mM DMPO as a spin trap (fig.1b).
- 3) hydroxyl radicals detection during the UV irradiation (λ=395 nm) of 10 mM H₂O₂ with 100 mM DMPO or 50 mM DEPMPO (fig 2).

Examples of experimental studies of oxidative processes by the EPR method using a benchtop CW X-band spectrometer demonstrate the broad possibilities of its practical applications. The method and equipment make it possible to effectively monitor the accumulation of free radicals in the sample of interest. In addition, due to the wide functionality of the SPINSCAN X EPR spectrometer, the EPR technique can also be used for deeper scientific research and study of the mechanisms of initiation and occurrence of redox reactions as well as in the creation of new oxidation-stable products.



Innovation
EPR spectrometer
SpinScan X

